MED-100

Báscula digital de cristal

Digital crystal scale

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V1.0

Manual de instrucciones

1015v

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Before using your new Steren's product, please read this manual to prevent any damage.

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The instructions of this manual are for reference about the product. There may be differences due to updates. Please check our web site (www.steren.com) to obtain the latest version of the instruction manual.

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- Use a soft and dry clothes to clean the scale.
- Don't use the scale if you have wet feet.
- Remove the battery if don't use it.
- Always place the scale onto a stable surface.
- Don't exceed the max charge capacity (150kg).
- Avoid dropping the unit as this may cause damage.

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Stay you in form, watching your weight with the new digital crystal scale. The scale is elegant, light and easy to use, you can adjust 3 types of units: Kg/Lb/St. Its platform is tempered crystal with 6 mm espesor, and have single sensors of high precision, as well as a screen LCD with 22mm numbers to easy reading, integrates the auto-off function, offering large life in their battery.

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

1. Remove the cover located in the rear part of the scale. Set the battery. Immediately the screen turn on.

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2. Select the weight unit (kg, lb or st) pressing the button as shown in the below image.

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1. Place feet on the scale and keeps as still as possible. The scale turns on immediately.

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2. Wait for a few seconds while the reading it's stabilized.

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3. The final reading appears flashing into the screen. Then you can retire from the scale.

ESPECIFICACIONES / SPECIFICATIONS

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Input: 3 V --- (one CR2016 battery)

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Maximum capacity: 330 lb

Graduación: 100 g

Increments for weight measurement: 0.2 lb

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Product: Digital crystal scale
Part number: MED-100
Brand: Steren

WARRANTY

This Steren product is warranted under normal usage against defects in workmanship and materials to the original purchaser for one year from the date of purchase.

CONDITIONS

- 1.- This warranty card with all the required information, invoice or purchase ticket, product box or package, and product, must be presented when warranty service is required.
- 2.- If the product stills on the warranty time, the company will repair it free of charge.
- 3.- The repairing time will not exceed 30 natural days, from the day the claim was received.
- 4.- Steren sell parts, components, consumables and accessories to customer, as well as warranty service, at any of the addresses mentioned later.

THIS WARRANTY IS VOID IN THE NEXT CASES:

If the product has been damaged by an accident, acts of God, mishandling, leaky batteries, failure to follow enclosed instructions, improper repair by unauthorized personnel, improper safe keeping, among others.

- a) The consumer can also claim the warranty service in the purchase establishment.
- b) If you lose the warranty card, we can reissue it, if you show the invoice or purchase ticket.

RETAILER INFORMATION

Name of the retailer _____
Address _____
Product _____
Brand _____
Serial number _____
Date of delivery _____

In case your product fails or have questions, please contact your nearest dealer. If you are in Mexico, please give a call to our Call Center.

01 800 500 9000

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ELECTRÓNICA STEREN DE MONTERREY, S.A.
ELECTRÓNICA STEREN DE TIJUANA, S.A. DE C.V.